


Product Name: Recombinant Human M6PR Protein
Catalog #: PEH2973



Summary

Name	MPRD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P20645
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TCDLVGEK GK ESEKELALVK RLKPLFNKSF ESTVGQGS DT YIYIFRVCRE AGNHTSGAGL VQINKSNGKE

image



Background

Alternative Names Cation-dependent mannose-6-phosphate receptor (CD Man-6-P receptor) (CD-MPR) (46 kDa mannose 6-phosphate receptor) (MPR 46)

Note

For Research Use Only , Not for Diagnostic Use