


Product Name: Recombinant Human ITGAL Protein
Catalog #: PEH2974



Summary

Name	ITAL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P20701
Host	E.coli
Species	Human
Molecular Mass	64 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AGRHFGYRVL QVGNGVIVGA PEGNSTGSL YQCQSGTGHC LPVTLRGSNY TSKYLGMTLA TDPTDGSILA CDPGLSRTCD QNTYLSGLCY LFRQNLQGPM LQGRPGFQEC IKGNVDLVFL FDGSMSLQPD EFQKILDFMK DVMKKLSNTS YQFAAVQFST SYKTEFDFSD YVKRKDPDAL LKHVKHMLLL TNTFGAINYV ATEVFREELG ARPDATKVL IITDGEATDS GNIDAAKDII RYIIGIGKHF QTKESQETLH KFASKPASEF VKILDTFEKL KDLFTELQKK IYVIEGTSKQ DLTSFNMELS SSGISADLSR GHAVVGAVGA KDWAGGFLDL KADLQDDTFI GNEPLTPEVR AGYLGYYVTW LPSRQKTSLL ASGAPRYQHM GRVLLFQEPQ GGGHWSQVQT IHGTQIGSYF GGELCGVDVD QDGETELLI GAPLFYGEQR GGRVFIYQRR QLGFEVSEL QGDPGYPLGR FGEAITALD INGDLVDVA VGAPLEEQGA VYIFNGRHGG LSPQPSQRIE GTQVLSGIQW FGRSIHG VKD LEGDGLADVA VGAESQMIVL

image



Product Name: Recombinant Human ITGAL Protein
Catalog #: PEH2974



Background

Alternative Names

Integrin alpha-L (CD11 antigen-like family member A) (Leukocyte adhesion glycoprotein LFA-1 alpha chain) (LFA-1A) (Leukocyte function-associated molecule 1 alpha chain) (CD antigen CD11a)

Note

For Research Use Only , Not for Diagnostic Use