



Product Name: Recombinant Human DDC Protein
Catalog #: PEH2976



Summary

Name	DDC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P20711
Host	E.coli
Species	Human
Molecular Mass	48 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KEMVDYMAN Y MEGIEGRQVY PDVEPGYLRP LIPAAAPQEP DTFEDIINDV EKIIMPGVTH WHSPYFFAYF PTASSYPAML ADMLCGAIGC IGFSWAASPA CTELETVMMD WLGMLELPK AFLNEKAGEG GGVIQGSASE ATLVALLAAR TKVIHRLQAA SPELTQAAIM EKLVA YSSDQ AHSSVERAGL IGGVKLKAIP SDGNFAMRAS ALQEALERDK AAGLIPFFMV ATLGTTTCCS FDNLLEV GPI CNKEDIWLHV DAAYAGSAFI CPEFRHLLNG VEFADSFNFN PHKWLLVNFD CSAMWVKKRT DLTGAFRLDP TYLKHS HQDS GLITDYRHWQ IPLGRRFRSL KMWFVFRMYG VKGLQAYIRK HVQLSHEFES LVRQDPRFEI CVEVILGLVC FRLKGSNKVN

image



Background



Product Name: Recombinant Human DDC Protein
Catalog #: PEH2976

Alternative Names

Aromatic-L-amino-acid decarboxylase (AADC) (EC 4.1.1.28) (DOPA decarboxylase) (DDC)

Note

For Research Use Only , Not for Diagnostic Use