


Product Name: Recombinant Human ITGA6 Protein
Catalog #: PEH2985



Summary

Name	ITA6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P23229
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ITASVEIQEP SSRRRVNSLP EVLPILNSDE PKTAHIDVHF LKEGCGDDNV CNSNLKLEYK FCTREGNQDK FSYLPIQKGV PELVLKDQKD IALEITVTNS PSNPRNPTKD

image



Background

Alternative Names	Integrin alpha-6 (CD49 antigen-like family member F) (VLA-6) (CD antigen CD49f) [Cleaved into: Integrin alpha-6 heavy chain; Integrin alpha-6 light chain; Processed integrin alpha-6 (Alpha6p)]
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ITGA6 Protein
Catalog #: PEH2985

