


Product Name: Recombinant Human CD34 Protein
Catalog #: PEH2993



Summary

Name	CD34_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P28906
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EATTNITETT VKFTSTSVIT SVYGNTNSSV QSQTSVISTV FTT PANVSTP ETTLKPSLSP GNVSDLSTTS TSLATSPTKP YTSSSPILSD IKAEIKCSGI REVKLTQGIC LEQNKTSSEA EFKKDRGEGL ARVLCGEEQA

image



Background

Alternative Names Hematopoietic progenitor cell antigen CD34 (CD antigen CD34)

Note

For Research Use Only , Not for Diagnostic Use