


Product Name: Recombinant Human NOS1 Protein
Catalog #: PEH2998



Summary

Name	NOS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P29475
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AEMKDMGIQV DRDLDGKSHK PLPLGVENDR VFNDLWGKGN VPVVLNNPYS EKEQPPTSGK QSPTKNGSPS KCPRFLKVKN

image



Background

Alternative Names	Nitric oxide synthase, brain (EC 1.14.13.39) (Constitutive NOS) (NC-NOS) (NOS type I) (Neuronal NOS) (N-NOS) (nNOS) (Peptidyl-cysteine S-nitrosylase NOS1) (bNOS)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use