


Product Name: Recombinant Human TIA1 Protein
Catalog #: PEH3004



Summary

Name	TIA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P31483
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VGNLSRDVTE ALILQLFSQI GPCKNCKMIM DTAGNDPYCF VEFHEHRHAA AALAAMNGRK IMGKEVKVNW ATTPSSQKKD TSSSTVVSTQ RSQDHFHV FV GDLSPEITTE DIKAAFAPFG RISDARVVKD MATGKSKGYG FVSFFNKWDA ENAIQQMGGQ WLGGRQIRTN WATRKPPAPK STYESNTKQL SYDEVVNQSS PSNCTVYCGG VTSGLTEQLM RQTFSPFGQI MEIRVFPDKG YSFVRFNSHE SAAHAIVSVN GTTIEGHVVK CYWGKETLDM INPVQQQNQI GYPQPYGQWG QWYGNAQQIG QYMPNGWQVP AYGMYGQAWN QQGFNQTQSS APWMGPNYGV QPPQGQNGSM LPNQPSGYRV

image



Background

Alternative Names Nucleolysin TIA-1 isoform p40 (RNA-binding protein TIA-1) (T-cell-restricted)

Product Name: Recombinant Human TIA1 Protein
Catalog #: PEH3004



intracellular antigen-1) (TIA-1) (p40-TIA-1)

Note

For Research Use Only , Not for Diagnostic Use