


Product Name: Recombinant Human HNRNPH3 Protein
Catalog #: PEH3005



Summary

Name	HNRH3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P31942
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NDASDGTVRL RGLPFGCSKE EIVQFFQGLE IVPNGITLTM DYQGRSTGEA FVQFASKEIA ENALGKHKER IGHRYIEIFR SSRSEIKGFY DPPRLLGQR PGPYDRPIGG RGGYYGAGRG SMYDRMRRGG DGYDGGYGGF DDYGGYNNYG YGNDGFDDRM

image



Background

Alternative Names Heterogeneous nuclear ribonucleoprotein H3 (hnRNP H3) (Heterogeneous nuclear ribonucleoprotein 2H9) (hnRNP 2H9)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human HNRNPH3 Protein
Catalog #: PEH3005

