


Product Name: Recombinant Human CYP2C19 Protein
Catalog #: PEH3009



Summary

Name	CP2CJ_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P33261
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HQESMDINNP RDFIDCFLIK MEKEKQNQQS EFTIENLVIT AADLLGAGTE TTSTTLRYAL LLLKHPEVT

image



Background

Alternative Names	Cytochrome P450 2C19 (EC 1.14.14.1) ((R)-limonene 6-monooxygenase) (EC 1.14.14.53) ((S)-limonene 6-monooxygenase) (EC 1.14.14.51) ((S)-limonene 7-monooxygenase) (EC 1.14.14.52) (CYP11C17) (CYP11C19) (Cytochrome P450-11A) (Cytochrome P450-254C) (Fenbendazole monooxygenase (4'-hydroxylating)) (EC 1.14.14.75) (Mephenytoin 4-hydroxylase)
--------------------------	---

Note

Product Name: Recombinant Human CYP2C19 Protein
Catalog #: PEH3009



For Research Use Only , Not for Diagnostic Use