



Product Name: Recombinant Human SDC2 Protein
Catalog #: PEH3011



Summary

Name	SDC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P34741
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDVESPELTT SRPLPKILLT SAAPKVETT LNIQNKIPAQ TKSPEETDKE KVHLSDSERK MDPAEEDTNV YTEKHSDSLF

image



Background

Alternative Names Syndecan-2 (SYND2) (Fibroglycan) (Heparan sulfate proteoglycan core protein) (HSPG) (CD antigen CD362)

Note

For Research Use Only , Not for Diagnostic Use