


Product Name: Recombinant Human AGL Protein
Catalog #: PEH3023



Summary

Name	GDE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P35573
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AVTRHPSIH QSVVAVSRTA FRNPKTSFYS KEVPQMCIPG KIEEVVLEAR TIERNTKPYR KDENSINGTP DITVEIREHI QLNESKIVKQ AGVATKGPNE

image



Background

Alternative Names	Glycogen debranching enzyme (Glycogen debrancher) [Includes: 4-alpha-glucanotransferase (EC 2.4.1.25) (Oligo-1,4-1,4-glucantransferase); Amylo-alpha-1,6-glucosidase (Amylo-1,6-glucosidase) (EC 3.2.1.33) (Dextrin 6-alpha-D-glucosidase)]
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human AGL Protein
Catalog #: PEH3023

