


Product Name: Recombinant Human SOX6 Protein
Catalog #: PEH3024



Summary

Name	SOX6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P35712
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TQDLTSREKE EGSDQHVASH LPLHPIMHMK PHSEELPTLV STIQQDADWD SVLSSQQRME SENNKLCSLY SFRNTSTSPH KPDEGSRDRE IMTSVTFGTP ERRKGLADV VDTLKQKKLE EMTRTEQEDS SCMEKLLSKD WKEKMERLNT SELLGEIKGT PESLAEKERQ LSTMITQLIS LREQLLAAHD EQKLAASQI EKQRQQMDLA RQQQEIQARQ QQQLLQQQHK INLLQQQIQV

image



Background

Alternative Names Transcription factor SOX-6

Note

Product Name: Recombinant Human SOX6 Protein
Catalog #: PEH3024



For Research Use Only , Not for Diagnostic Use