


Product Name: Recombinant Human SOX18 Protein
Catalog #: PEH3025



Summary

Name	SOX18_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P35713
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AFMVWAKDER KRLAQQNPD L HNAVL SKMLG KAWKELNAAE KRPFVEEAER LRVQH LRDHP NYKYRPRRKK QARKARRLEP

image



Background

Alternative Names Transcription factor SOX-18

Note

For Research Use Only , Not for Diagnostic Use