


Product Name: Recombinant Human GRIK1 Protein
Catalog #: PEH3035



Summary

Name	GRIK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P39086
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSIEKWSME RLQAPRPET GLLDGMMTTE AALMYDAVYM VAIASHRASQ LTVSSLQCHR HKPWRLGPRF MNLKEARWD GLTGHTFNK TNGLRKDFDL DIISLKEEGT EKAAGEVSKH LYKVWKKIGI WNSNSGLNMT DSNKDKSSNI TDSLANTLI VTTILEEPYV MYRKSDKPLY GNDRFEGYCL DLLKELSNIL GFIYDVKLVP DGKYGAQNDK GEWNGMVKEL IDHRADLAVA

image



Background

Alternative Names Glutamate receptor ionotropic, kainate 1 (GluK1) (Excitatory amino acid receptor 3) (EAA3) (Glutamate receptor 5) (GluR-5) (GluR5)

Note

Product Name: Recombinant Human GRIK1 Protein
Catalog #: PEH3035



For Research Use Only , Not for Diagnostic Use