

Product Name: Recombinant Human IL6ST Protein
Catalog #: PEH3036



Summary

| | |
|------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Name | IL6RB_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | P40189 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 32 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | IWPNVPDPSK SHIAQWSPHT PPRHNFNSKD QMYSDGNFTD VSVVEIEAND KKFPEDLKS LDLFKKEKIN TEGHSSGIGG SSCMSSSRPS ISSSDENESS QNTSSTVQYS TVVHSGYRHQ VPSVQVFSRS ESTQPLLDSE ERPEDLQLVD HVDGGDGILP RQQYFKQNCS QHESSPDISH FERSKQVSSV NEEDFVRLKQ QISDHISQSC GSGQMCMFQE VSAADAFGPG TEGQVERFET VGMEAATDEG MPKSYLPQTV |

image



Background

| | |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Alternative Names | Interleukin-6 receptor subunit beta (IL-6 receptor subunit beta) (IL-6R subunit beta) (IL-6R-beta) (IL-6RB) (CDw130) (Interleukin-6 signal transducer) (Membrane glycoprotein 130) (gp130) (Oncostatin-M receptor) |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

Product Name: Recombinant Human IL6ST Protein
Catalog #: PEH3036



subunit alpha) (CD antigen CD130)

Note

For Research Use Only , Not for Diagnostic Use