


Product Name: Recombinant Human PSMB10 Protein
Catalog #: PEH3041



Summary

Name	PSB10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P40306
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TTIAGLVFQD GVLGADTRA TNSVVDKSK CEKIHFIAPK IYCCGAGVAA DAEMTTRMVA SKMELHALST GREPRVATVT RILRQTLFRY QGHV GASLIV GGVDLTGPQL YGVHPHGSYS RLPFTALGSG QDAALAVLED RFQPNMTLEA AQGLLVEAVT AGILGDLGSG GNVDACVITK TGAKLLRRTLSP SPTPEVKRSG RYHFVPGTTA VLTQTVKPLT LELVEETVQA MEVE

image



Background

Alternative Names	Proteasome subunit beta type-10 (EC 3.4.25.1) (Low molecular mass protein 10) (Macropain subunit MECl-1) (Multicatalytic endopeptidase complex subunit MECl-1) (Proteasome MECl-1) (Proteasome subunit beta-2i)
--------------------------	---

Product Name: Recombinant Human PSMB10 Protein
Catalog #: PEH3041



Note

For Research Use Only , Not for Diagnostic Use