


Product Name: Recombinant Human CASR Protein
Catalog #: PEH3043



Summary

Name	CASR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P41180
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EVRCSAAHA FKVAARATLR RSNVSRKRSS SLGGSTGSTP SSSISSKSNS EDPFPQPERQ KQQQPLALTQ QEQQQQPLTL PQQQRSQQQP RCKQKVIFGS GTVTFSLSFD EPQKNAMHR NSTHQNSLEA QKSSDTLTRH EPLLPLQCGE TDLTLTVQET GLQGPVGGDQ RPEVEDPEEL SPALVVSSSQ SFVISGGGST

image



Background

Alternative Names	Extracellular calcium-sensing receptor (CaR) (CaSR) (hCasR) (Parathyroid cell calcium-sensing receptor 1) (PCaR1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CASR Protein
Catalog #: PEH3043

