


Product Name: Recombinant Human GRIA3 Protein
Catalog #: PEH3051



Summary

Name	GRIA3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P42263
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NEYIEQRKPC DTMKVGGNLD SKGYGVATPK GSALRNAVNL AVLKLNEQGL LDKLNKWWY DKGECGSGGG

image



Background

Alternative Names Glutamate receptor 3 (GluR-3) (AMPA-selective glutamate receptor 3) (GluR-C) (GluR-K3) (Glutamate receptor ionotropic, AMPA 3) (GluA3)

Note

For Research Use Only , Not for Diagnostic Use