

Product Name: Recombinant Human DCC Protein
Catalog #: PEH3053



Summary

Name	DCC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P43146
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKRATHSAGK RKGSKDLRP PDLWIHHEEM EMKNIEKPSG TDPAGRDSPI QSCQDLTPVS HSQSETQLGS

image



Background

Alternative Names Netrin receptor DCC (Colorectal cancer suppressor) (Immunoglobulin superfamily DCC subclass member 1) (Tumor suppressor protein DCC)

Note

For Research Use Only , Not for Diagnostic Use