



Summary

Name	KCNJ4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P48050
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGEYLPLDQR DLNMGYDIGL DRIFLVSPII IVHEIDEDSP LYGMGKEELE SEDFEIVVIL EGMVEATAMT TQARSSYLAS

image



Background

Alternative Names	Inward rectifier potassium channel 4 (HIRK2) (HRK1) (Hippocampal inward rectifier) (HIR) (Inward rectifier K(+)) channel Kir2.3) (IRK-3) (Potassium channel, inwardly rectifying subfamily J member 4)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use