



Product Name: Recombinant Human SOX9 Protein
Catalog #: PEH3064



Summary

Name	SOX9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P48436
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NKPHVKRPMN AFMVWAQAAR RKLADQYPHL HNAELSKTLG KLRLLNESE KRPFVEEAER LRVQHKKDHP DYKYQPRRRK SVKNGQAEAE EATEQTHISP NAIFKALQAD SPHSSSGMSE VHSPGEHSGQ SQGPPTPPTT PKTDVQPGKA DLKREGRPLP EGGRQPPIDF

image



Background

Alternative Names Transcription factor SOX-9

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SOX9 Protein
Catalog #: PEH3064

