


Product Name: Recombinant Human SLC4A3 Protein
Catalog #: PEH3066



Summary

Name	B3A3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P48751
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HSSSSPSRA RASRLAGEKS RPWSPSASYD LRERLCPGSA LGNPGGPEQQ VPTDEAEAQM LGSADLDDMK SHRLEDNPGV RRHLVKKPSR TQGGRGSPSG LAPILRRKKK KKKLDRRPHE VFVELNELML DRSQEPHWRE

image



Background

Alternative Names	Anion exchange protein 3 (AE 3) (Anion exchanger 3) (CAE3/BAE3) (Cardiac/brain band 3-like protein) (Neuronal band 3-like protein) (Solute carrier family 4 member 3)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SLC4A3 Protein
Catalog #: PEH3066

