


Product Name: Recombinant Human TRPC1 Protein
Catalog #: PEH3067



Summary

Name	TRPC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P48995
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VMALKDVREV KEENTLNEKL FLLACDKGDY YMVKKILEEN SSGDLNINCV DVLGRNAVTI TIENENLDIL

image



Background

Alternative Names	Short transient receptor potential channel 1 (TrpC1) (Transient receptor protein 1) (TRP-1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use