


Product Name: Recombinant Human METAP2 Protein
Catalog #: PEH3071



Summary

Name	MAP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P50579
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AGEQEPDKES GASVDEVARQ LERSALEDKE RDEDDGDG DGDGATGKKK KKKKKKRGPK VQTDPPSVPI CDLYPNGVFP KGQCEYPPT QDGRTAAWRT TSEKKALDQ ASEEIWNDFR EAAEAHRQVR KYVMSWIKPG

image



Background

Alternative Names	Methionine aminopeptidase 2 (MAP 2) (MetAP 2) (EC 3.4.11.18) (Initiation factor 2-associated 67 kDa glycoprotein) (p67) (p67eIF2) (Peptidase M)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use