


Product Name: Recombinant Human FXR1 Protein
Catalog #: PEH3073



Summary

| | |
|------------------------------|--|
| Name | FXR1_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | P51114 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 31 kDa |
| Buffer | Lyophilized from Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | EEFVVREDLM GLAIGTHGSN IQQARKVPGV TAIELDEDTG TFRIYGESAD AVKKARGFLE FVEDFIQVPR NLVGKVIKGN GKVIQEIVDK SGVVRVRIEG DNENKLPRED GMVPFVFGVT KESIGNVQVL LEYHIAYLKE VEQLRMERLQ IDEQLRQIGS RSYSGRGRGR RGPNYTSGYG TNSELSNPSE TESERKDELS DWSLAGEDDR DSRHQRDSRR RPGGRGRSVS GGRGRGGPRG GKSSISSVLK DPDSNPYSLL |

image



Background

Alternative Names Fragile X mental retardation syndrome-related protein 1 (hFXR1p)

Product Name: Recombinant Human FXR1 Protein
Catalog #: PEH3073



Note

For Research Use Only , Not for Diagnostic Use