



Product Name: Recombinant Human GFRA1 Protein
Catalog #: PEH3084



Summary

Name	GFRA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P56159
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QCVAGKETNF SLASGLEAKD ECRSAMEALK QKSLYNCRCK RGMKKEKNCL RIYWSMYQSL QGNDLLEDSP YEPVNSRLSD IFRVVPFISD VFQQVEHIPK GNNCLDAAKA CNLDDICKKY RSAYITPCTT SVSNDVCNRR KCHKALRQFF DKVPAKHSYG MLFCSCRDIA CTERRRQTIV PVCSYEEREK PNCLNLQDSC KTNYICRSRL ADFFTNCQPE SRSVSSCLKE

image



Background

Alternative Names	GNDF family receptor alpha-1 (GNDF receptor alpha-1) (GDNFR-alpha-1) (GFR-alpha-1) (RET ligand 1) (TGF-beta-related neurotrophic factor receptor 1)
--------------------------	---

Product Name: Recombinant Human GFRA1 Protein
Catalog #: PEH3084



Note

For Research Use Only , Not for Diagnostic Use