

Product Name: Recombinant Human TPTE Protein
Catalog #: PEH3085



Summary

Name	TPTE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P56180
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AGVIIELGPN DSPQTSEFKG ATEEAPAKES PHTSEFKGAA RVSPIESVL ARLSKFEVED AENVASYDSK

image



Background

Alternative Names	Putative tyrosine-protein phosphatase TPTE (EC 3.1.3.48) (Cancer/testis antigen 44) (CT44) (Transmembrane phosphatase with tensin homology) (Tumor antigen BJ-HCC-5)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use