


Product Name: Recombinant Human MAZ Protein
Catalog #: PEH3086



Summary

Name	MAZ_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P56270
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VCQQRFKRKD RMSYHVRSHD GAVHKPYNCS HCGKSFSRPD HLNSHVRQVH STERPFKCEK CEAAFATKDR LRAHTVRHEE KVPCHVCGKM LSSAYISDHM KVHSQGP HHV CELCNKGTGE VCPMAAAAAA

image



Background

Alternative Names Myc-associated zinc finger protein (MAZI) (Pur-1) (Purine-binding transcription factor) (Serum amyloid A-activating factor-1) (SAF-1) (Transcription factor Zif87) (ZF87) (Zinc finger protein 801)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MAZ Protein
Catalog #: PEH3086

