



Product Name: Recombinant Human HDAC4 Protein
Catalog #: PEH3089



Summary

Name	HDAC4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P56524
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AEFQRQHEQL SRQHEAQLHE HIKQQQEMLA MKHQQELLEH QRKLERHRQE QELEKQHREQ KLQQLKNKEK GKESAVASTE VKMKLQEFVL NKKKALAHRN LNHCISDPR YWYGKTQHSS LDQSSPPQSG VSTSYNHPVL GMYDAKDDFP LRKTASEPNL KLSRSLKQKV AERRSSPLL RKGDPVV TAL KKRPLDVTDS ACSSAPGSGP SSPNNSGGSV SAENGIAPAV

image



Background

Alternative Names Histone deacetylase 4 (HD4) (EC 3.5.1.98)

Note

Product Name: Recombinant Human HDAC4 Protein
Catalog #: PEH3089



For Research Use Only , Not for Diagnostic Use