


Summary

Name	SRP54_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P61011
Host	E.coli
Species	Human
Molecular Mass	49 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TSALRSLSNA TIINEEVLNA MLKEVCTALL EADVNIKLVK QLRENVKSAI DLEEMASGLN KRKMIQHAFV KELVKLVDPG VKAWTPTKGG QNVIMFVGLQ GSGKTTTCSK LAYYYQRKGW KTCLICADTF RAGAFDQLKQ NATKARIPFY GSYTEMDPVI IASEGVEKFK NENFEIIIIVD TSGRHKQEDS LFEEMLQVAN AIQPDNIVYV MDASIGQACE AQAKAFKDKV DVASVIVTKL DGHAKGGGAL SAVAATKSPI IFIGTGEHID DFEPFKTQPF ISKLLGMGDI EGLIDKVNEL KLDDNEALIE KCLKHGQFTLR DMYEQFQNM KMGPFSQLG MIPGFGTDFM SKGNEQESMA RLKCLMTIMD SMNDQELDST DGAKVFSKQP GRIQRVARG GVSTRDVQEL LTQYTKFAQM

image



Background

Product Name: Recombinant Human SRP54 Protein
Catalog #: PEH3098



Alternative Names

Signal recognition particle 54 kDa protein (SRP54) (EC 3.6.5.-)

Note

For Research Use Only , Not for Diagnostic Use