


Product Name: Recombinant Human LMO4 Protein
Catalog #: PEH3099



Summary

Name	LMO4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P61968
Host	E.coli
Species	Human
Molecular Mass	57 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KRCAGCGGKI ADRFLLYAMD SYWHSRCLKC SCCQAQLGDI GTSCYTKSGM ILCRNDYIRL FGNSGACSAC GQSIPASELV MRAQGNVYHL KCFTCSTCRN RLVPGDRFHY INGS LFCEHD RPTALINGHL

image



Background

Alternative Names	LIM domain transcription factor LMO4 (Breast tumor autoantigen) (LIM domain only protein 4) (LMO-4)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use