


Product Name: Recombinant Human SNRPE Protein
Catalog #: PEH3103



Summary

Name	RUXE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P62304
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAYRGQGQKV QKVMVQPINL IFRYLQNRSR IQVWLYEQVN MRIEGCIIGF DEYMNVLVDD AEEIHSKTKS RKQLGRIMLK GDNITLLQSV SN

image



Background

Alternative Names Small nuclear ribonucleoprotein E (snRNP-E) (Sm protein E) (Sm-E) (SmE)

Note

For Research Use Only , Not for Diagnostic Use