


Product Name: Recombinant Human SNRPD3 Protein
Catalog #: PEH3107



Summary

Name	SMD3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P62318
Host	E.coli
Species	Human
Molecular Mass	54 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HEAEGHIVTC ETNTGEVYRG KLIEAEDNMN CQMSNITVTY RDGRVAQLEQ VYIRGSKIRF LILPDMLKNA PMLKSMKNKN QGSGAGRGKA

image



Background

Alternative Names Small nuclear ribonucleoprotein Sm D3 (Sm-D3) (snRNP core protein D3)

Note

For Research Use Only , Not for Diagnostic Use