


**Product Name: Recombinant Human OLR1 Protein**  
**Catalog #: PEH3115**



## Summary

<b>Name</b>	OLR1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P78380
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	12 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	SQVSDLLTQE QANLTHQKKK LEGQISARQQ AEEASQSEEN ELKEMIETLA RKLNEKSKEQ MELHHQNLNL QETLKRIVANC SAPCPQDWIW

## image



## Background

<b>Alternative Names</b>	Oxidized low-density lipoprotein receptor 1 (Ox-LDL receptor 1) (C-type lectin domain family 8 member A) (Lectin-like oxidized LDL receptor 1) (LOX-1) (Lectin-like oxLDL receptor 1) (hLOX-1) (Lectin-type oxidized LDL receptor 1) [Cleaved into: Oxidized low-density lipoprotein receptor 1, soluble form]
--------------------------	--

## Note

**Product Name: Recombinant Human OLR1 Protein**  
**Catalog #: PEH3115**



---

For Research Use Only , Not for Diagnostic Use