



Product Name: Recombinant Human NUCB2 Protein
Catalog #: PEH3118



Summary

Name	NUCB2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P80303
Host	E.coli
Species	Human
Molecular Mass	48 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KTKVQNIHPV ESAKIEPPDT GLYYDEYLKQ VIDVLETDKH FREKLQKADI EEIKSGRLSK ELDLVSHHVR TKLDELKRQE VGRLRMLIKA KLDSLQDIGM DHQALLKQFD HLNHLNPDKF ESTDLDMLIK AATSDLEHYD KTRHEEFKKY EMMKEHERRE YLKTLENEKR KEEESKFEEM KKKHENHPKV NHPGSKDQLK EVWEETDGLD PNDFDPKTFE KLHDVNSDGF LDEQELEALF TKELEKVYDP KNEEDDMVEM EEERLRMREH VMNEVDTNKD RLVTLLEFLK ATEKKEFLEP DSWETLDQQQ FFTEELKEY ENIIALQENE LKKKADELQK QKEELQRQHD QLEAQKLEYH QVIQQMEQKK LQQGIPPSGP

image



Background

Alternative Names	Nucleobindin-2 (DNA-binding protein NEFA) (Epididymis secretory protein Li
--------------------------	--

Product Name: Recombinant Human NUCB2 Protein
Catalog #: PEH3118



109) (Gastric cancer antigen Zg4) (Prepronesfatin) [Cleaved into: Nesfatin-1]

Note

For Research Use Only , Not for Diagnostic Use