


Product Name: Recombinant Human MPP1 Protein
Catalog #: PEH3121



Summary

Name	EM55_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q00013
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GGSMHTALSD LYLEHLLQKR SRPEAVSHPL NTVTEDMYTN GSPAGSPAQ VKGQEVKVR LIQFEKVTEE PMGITLKLNE KQSCTVARIL HGGMIHRQGS LHVGDEILEI NGTNVTNHSV DQLQKAMKET KGMISLKVIP NQQSRLPALQ MFMRAQFDYD PKKDNLIPCK EAGLK FATGD IIQIINKDDS NWWQGRVEGS SKESAGLIPS PELQEW RVAS MAQSAPSEAP SCSPFGKKKK YKDKYLAKHS SIFDQLDVVS YEEVVRLPAF

image



Background

Alternative Names	55 kDa erythrocyte membrane protein (p55) (Membrane protein, palmitoylated 1)
--------------------------	---

Product Name: Recombinant Human MPP1 Protein
Catalog #: PEH3121



Note

For Research Use Only , Not for Diagnostic Use