

Product Name: Recombinant Human CDR2 Protein
Catalog #: PEH3126



Summary

Name	CDR2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q01850
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DEEENEHLKK TVTMLQAQLS LERQKRV TME EEYGLVLKEN SELEQQLGAT GAYRARALEL EAEVAEMRQM LQSEHPFVNG VEKLV PDSLY VPFKEPSQSL LEEMFLT VPE SHRKPLKRSS SETILSSLAG SDIVKGHEET CIRRAKAVKQ RGISLLHEVD TQYSALKVKY EELLKKCQEE QDSL SHKAVQ TSRAAAKDLT GVNAQSEPVA SGWELASVNP EPVSSPTTPP EYKALFKEIF SCIKKTKQEI DEQRTKYRSL

image



Background

Alternative Names	Cerebellar degeneration-related protein 2 (Major Yo paraneoplastic antigen) (Paraneoplastic cerebellar degeneration-associated antigen)
--------------------------	--

Product Name: Recombinant Human CDR2 Protein
Catalog #: PEH3126



Note

For Research Use Only , Not for Diagnostic Use