


Product Name: Recombinant Human CNTN2 Protein
Catalog #: PEH3128



Summary

Name	CNTN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q02246
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PGGVVVRDIG DTTIQLSWSR GFDNHSPIAK YTLQARTPPA GKWKQVRTNP ANIEGNAETA QVLGLTPWMD YEFRVIASNI LGTGEPSGPS SKIRTREAAP SVAPSGLSGG GGAPGELIVN WTPMSREYQN GDGFGYLLSF RRQGSTHWQT ARVPGADAQY FVYSNESVRP YTPFEVKIRS YNRRGDGPES LTALVYSAEE

image



Background

Alternative Names	Contactin-2 (Axonal glycoprotein TAG-1) (Axonin-1) (Transient axonal glycoprotein 1) (TAX-1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CNTN2 Protein
Catalog #: PEH3128

