



Product Name: Recombinant Human DST Protein
Catalog #: PEH3132



Summary

Name	DYST_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q03001
Host	E.coli
Species	Human
Molecular Mass	67 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGFSIQEKYV AADTLYSQIK EDVKKRAVAL DEAISQSTQF HDKIDQILES LERIVERLRQ PPSISAEVEK IKEQISENKN VSVDMEKLQP LYETLKQRGE EMIARSGGTD KDISAKAVQD KLDQMVFIWE NIHTLVEERE AKLLDVMELA EKFWCDHMSL IVTIKDTQDF IRDLEDPGID PSVVKQQQEA AETIREEIDG LQEELDIVIN LGSELIAACG EPDKPIVKKS IDELSAWDS LNKAWKDRID KLEEAMQAAV QYQDGLQAVF DWVDIAGGKL ASMSPIGTDL ETVKQQIEEL KQFKSEAYQQ QIEMERLNHQ AELLKKVTE ESKHTVQDP LMELKLIWDS LEERIIINRQH KLEGALLALG QFQHALDELL AWLTHTEGLL SEQKPVGGDP KAIEIELAKH HVLQNDVLAH QSTVEAVNKA GNDLISSAG EEASNLQNKL EVLNQRWQNV LEKTEQRKQQ LDGALRQAKG FHGEIEDLQQ WLTDTERHLL ASKPLGGLPE TAKEQLNVHM EVCAAFEAKE ETYKSLMQKG QQMLARCPKS AETNIDQDIN NLKEKWESVE

image



Product Name: Recombinant Human DST Protein
Catalog #: PEH3132



Background

Alternative Names

Dystonin (230 kDa bullous pemphigoid antigen) (230/240 kDa bullous pemphigoid antigen) (Bullous pemphigoid antigen 1) (BPA) (Bullous pemphigoid antigen) (Dystonia musculorum protein) (Hemidesmosomal plaque protein)

Note

For Research Use Only , Not for Diagnostic Use