


Product Name: Recombinant Human CENPC Protein
Catalog #: PEH3133



Summary

Name	CENPC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q03188
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EFYLSVGSPS VLLDAKTSVS QNVIPSSAQK RETYTFENSV NMLPSSTEVS VKTKKRLNFD DKVMLKKIEI DNKVSDEEDK

image



Background

Alternative Names	Centromere protein C (CENP-C) (Centromere autoantigen C) (Centromere protein C 1) (CENP-C 1) (Interphase centromere complex protein 7)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use