


Product Name: Recombinant Human ICA1 Protein
Catalog #: PEH3136



Summary

Name	ICA69_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q05084
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PYEFTTLKSL QDPMKKLVEK EEKKKINQQE STDAAVQEPS QLISLEENQ RKESSSFKTE DGKSILSALD KGSTHTACSG PIDELDMKS

image



Background

Alternative Names	Islet cell autoantigen 1 (69 kDa islet cell autoantigen) (ICA69) (Islet cell autoantigen p69) (ICAp69) (p69)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use