


Product Name: Recombinant Human CKAP4 Protein
Catalog #: PEH3142



Summary

Name	CKAP4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q07065
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HQDFSRQREE LGQGLQGVEQ KVQSLQATFG TFESILRSSQ HKQDLTEKAV KQGESEVSRI SEVLQKLQNE ILKDLSDGIH VVKDARERDF TSLENTVEER LTELTKSIND NIAIFTEVQK RSQKEINDMK AKVASLEESE GNKQDLKALK EAVKEIQTSA KSREWDMEAL RSTLQTMESD IYTEVRELVS LKQEQQAFKE AADTERLALQ ALTEKLLRSE ESVSRLPEEI RRLEEELRQL KSDSHGPKED GGFRHSEAFE ALQQKSQGLD SRLQHVEDGV LSMQVASARQ TESLESLLSK SQEHEQRLAA LQGRLEGLGS

image



Background

Alternative Names Cytoskeleton-associated protein 4 (63-kDa cytoskeleton-linking membrane protein) (Climp-63) (p63)

Product Name: Recombinant Human CKAP4 Protein
Catalog #: PEH3142



Note

For Research Use Only , Not for Diagnostic Use