


Product Name: Recombinant Human LRP1 Protein
Catalog #: PEH3143



Summary

Name	LRP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q07954
Host	E.coli
Species	Human
Molecular Mass	63 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HGLAVDWVSR NLFWTSYDTN KKQINVARLD GSFKNAVVGQ LEQPHGLVVH PLRGKLYWTD GDNISMANMD GSNRTLLFSG QKGPVGLAID FPESKLYWIS SGNHTINRCN LDGSGLEVID AMRSQLGKAT ALAIMGDKLW WADQVSEKMG TCSKADGSGS VVLRNSTTLV MHMKVYDESI QLDHKGTNPC SVNNGDCSQL CLPTSETTRS CMCTAGYSLR SGQQACEGVG SFLLYSVHEG IRGIPLDPND KSDALVPVSG TSLAVGIDFH AENDTIYWVD MGLSTISRAK RDQTWREDEVV TNGIGRVEGI AVDWIAGNIY WTDQGFVIE VARLNGSFRY VVISQGLDKP RAITVHPEKG YLFWTEWGQY PRIERSRLDG TERVVLVNVS ISWPNGISVD YQDGKLYWCD ARTDKIERID LETGENREV LSSNNMDMFS VSVFEDFIYW SDRTHANGSI KRGSKNATD SVPLRTGIGV QLKDIKFVNR DRQKGTNVCA VANGGCQQLC LYRGRGQRAC ACAHGMLAED GASCREYAGY LLYSERTILK

image



Background

Product Name: Recombinant Human LRP1 Protein
Catalog #: PEH3143



Alternative Names

Prolow-density lipoprotein receptor-related protein 1 (LRP-1) (Alpha-2-macroglobulin receptor) (A2MR) (Apolipoprotein E receptor) (APOER) (CD antigen CD91) [Cleaved into: Low-density lipoprotein receptor-related protein 1 85 kDa subunit (LRP-85); Low-density lipoprotein receptor-related protein 1 515 kDa subunit (LRP-515); Low-density lipoprotein receptor-related protein 1 intracellular domain (LRPICD)]

Note

For Research Use Only , Not for Diagnostic Use