


Product Name: Recombinant Human CNTN1 Protein
Catalog #: PEH3150



Summary

Name	CNTN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q12860
Host	E.coli
Species	Human
Molecular Mass	53 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PADIVVQFKD VYALMGQNVV LECFALGNPV PDIRWRKVLV PMPSTAEIST SGAVLKIFNI QLEDEGIYEC EAENIRGKDK HQARIYVQAF PEWVEHINDT EVDIGSDLYW PCVATGKPIP TIRWLKNGYA YHKGELRLYD VTFENAGMYQ CIAENTYGAI YANAELKILA LAPTFEMNPM KKKILAAKGG RVIIIECKPKA APKPKFSWSK GTEWLVN SSR ILIWEDGSLE INNITRNDGG IYTCFAENNR GKANSTGTLV ITDPTRIILA PINADITVGE NATMQCAASF DPALDLTFVW SFNGYVIDFN KENIH YQRNF MLDSNGELLI RNAQLKHAGR YTCTAQTIVD NSSASADLVV

image



Background

Alternative Names Contactin-1 (Glycoprotein gp135) (Neural cell surface protein F3)

Product Name: Recombinant Human CNTN1 Protein
Catalog #: PEH3150



Note

For Research Use Only , Not for Diagnostic Use