


Product Name: Recombinant Human GRIN2A Protein
Catalog #: PEH3153



Summary

Name	NMDE1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q12879
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QKDNLNYYVF QGQHPLTLNE SNPNTVEVAV STESKANSRP RQLWKKSVDS IRQDSLSQNP VSQRDEATAE NRTHSLKSPR

image



Background

Alternative Names	Glutamate receptor ionotropic, NMDA 2A (GluN2A) (Glutamate [NMDA] receptor subunit epsilon-1) (N-methyl D-aspartate receptor subtype 2A) (NMDAR2A) (NR2A) (hNR2A)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use