


Product Name: Recombinant Human ELAVL2 Protein
Catalog #: PEH3157



Summary

Name	ELAV2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q12926
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	METQLSNGPT CNNTANGPTT INNNCSSPVD SGNTEDSKTN LIVNYLPQNM TQEELKSLFG SIGEIESCKL VRDKITGQSL GYGFVNYIDP KDAEKAINTL NGLRLQTKTI KVSYPARSSA SIRDANLYVS GLPKTMTQKE LEQLFSQYGR IITSRILVDQ VTGISRGVGF IRFDKRIEAE EAIKGLNGQK PPGATEPITV KFANNPSQKT NQAILSQLYQ SPNRRYPGPL AQQAQRFRLD NLLNMAYGVK RFSPMTIDGM TSLAGINIPG HPGTGWCIFV YNLAPDADES ILWQMFGPFG AVTNVKVIRD FNTNKCKGFG FVTMTNYDEA AMAIASLNGY RLGDRVLLQVS FKTNKTHKA

image



Background

Alternative Names	ELAV-like protein 2 (ELAV-like neuronal protein 1) (Hu-antigen B) (HuB)
--------------------------	---

Product Name: Recombinant Human ELAVL2 Protein
Catalog #: PEH3157



(Nervous system-specific RNA-binding protein Hel-N1)

Note

For Research Use Only , Not for Diagnostic Use