


Product Name: Recombinant Human GPS1 Protein
Catalog #: PEH3159



Summary

Name	CSN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13098
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YEEIHRKLSE ATRSSLRELQ NAPDAIPESG VEPPALDTAW VEATRKALL KLEKLDLDL NYKGNSIKES IRRGHDDLGD

image



Background

Alternative Names	COP9 signalosome complex subunit 1 (SGN1) (Signalosome subunit 1) (G protein pathway suppressor 1) (GPS-1) (JAB1-containing signalosome subunit 1) (Protein MFH)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use