


Product Name: Recombinant Human PRDX4 Protein
Catalog #: PEH3161



Summary

Name	PRDX4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13162
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DHSLHLSKAK ISKPAPYWEG TAVIDGEFKE LKLTDIRGKY LVFFFYPLDF TFVCPTEIIA FGDRLEEFVS INTEVVACSV DSQFTHLAWI NTPRRQGGLG PIRIPLSDL THQISKDYGV YLEDSGHTLR GLFIIDDKGI LRQITLNDLP VGRSVDETLP LVQAFQYTDK HGEVCPAGWK PGSETIIPDP

image



Background

Alternative Names	Peroxiredoxin-4 (EC 1.11.1.24) (Antioxidant enzyme AOE372) (AOE37-2) (Peroxiredoxin IV) (Prx-IV) (Thioredoxin peroxidase A0372) (Thioredoxin-dependent peroxide reductase A0372) (Thioredoxin-dependent peroxiredoxin 4)
--------------------------	--

Product Name: Recombinant Human PRDX4 Protein
Catalog #: PEH3161



Note

For Research Use Only , Not for Diagnostic Use