


Product Name: Recombinant Human MMRN1 Protein
Catalog #: PEH3163



Summary

Name	MMRN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13201
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LSTEQVSDQK NAPAASVSN NVTEYMSTLH ENIKKQSLMM LQMFEDLHIQ ESKINNLTVS LEMEKESLRG ECEDMLSKCR NDFKFQLKDT EENLHVLNQT LAEVLFPM DN KMDKMSEQLN DLT YDMEILQ PLLEQGASLR QTMTYEQPKE AIVIRKKIEN LTS AVNSLNF IIKELTKRHN LLRNEVQGRD DALERRINEY ALEMEDGLNK TMTIINNAID FIQDN YALKE TLSTIKDNSE IHHKCTSDME

image



Background

Alternative Names	Multimerin-1 (EMILIN-4) (Elastin microfibril interface located protein 4) (Elastin microfibril interfacier 4) (Endothelial cell multimerin) [Cleaved into: Platelet glycoprotein Ia*; 155 kDa platelet multimerin (p-155) (p155)]
--------------------------	---

Product Name: Recombinant Human MMRN1 Protein
Catalog #: PEH3163



Note

For Research Use Only , Not for Diagnostic Use