


Product Name: Recombinant Human XRCC4 Protein
Catalog #: PEH3171



Summary

Name	XRCC4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13426
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VSEPSITHFL QVSWEKTLES GFVITLTDGH SAWTGTVSES EISQEADDMA MEKGKYVGEL RKALLSGAGP ADVYTFNFSK ESCYFFFEKN LKDVSFRLGS FNLEKVENPA EVIRELICYC LDTIAENQAK NEHLQKENER LLRDWNDVQG RFEKCVSAKE ALETDLYKRF ILVLNEKGTK IRSLHNKLLN AAQEREKDIK QEGETAICSE MTADRPVYD ESTDEESENQ TDLGLASAA VSKDDSISS LDVTDIAPSR KRRQRMQRNL GTEPKMAPQE NQLQEKENS R PDSSLPETSK KEHISAENMS LETLRNSSPE

image



Background

Alternative Names DNA repair protein XRCC4 (X-ray repair cross-complementing protein 4)

Product Name: Recombinant Human XRCC4 Protein
Catalog #: PEH3171



Note

For Research Use Only , Not for Diagnostic Use